




simple&smart

KISS 

Keep It Simple & Smart



SIMPLE inside
less assistance need
easy installation
LESS electronics
for a **MORE RELIABLE UNIT**
That's why we are
simple&smart




1 CHARACTERISTICS



CHARACTERISTICS

COMFY & COZY

improving ergonomics



The **SEAT** and the **BACKREST** have been designed to **improve the patient's comfort and stability** on the chair. The particular ergonomics allows the doctor to work in a more suitable position. Working in 3 o'clock and 9 o'clock position will be easier.

The structure of the **BACKREST** is in metal and its extra thin design improves the ergonomics of the seat and makes the ap-


proach of the doctor to the patient much more easier, in this way the doctor can reach a **better working position**.

The design of the **CHAIR BASE** grants a wider movement of the seat. It is thin and allows the chair to reach a **lower position**, to grant a wider working space and improve the doctor's comfort during the daily practice.

CHARACTERISTICS

GOOD TO THE BONE

resistent structure, performing features



The **STRUCTURE** of the unit is built to last in time, solid and resistant grant a perfect stability and durability.

The **WATER UNIT** can rotate to allow the doctor to work at 3 o'clock or to allow a better positioning of the assistant. The water unit's internal structure and the most fragile parts of the dental unit are made in **aluminum**. This feature allows to reduce the chair weight, making the movement easier to manage.

The SS_KISS is equipped with a **DISINFECTION SYSTEM**. There is a tank inside

the water unit for the disinfection of the water lines that can be filled with a proper product and activated by a special switch located inside the water unit, to avoid any possible accidental activation during the daily work.

There are two separated **AIRTIGHT FILTERS**, one for each suction cannula.

The SS_KISS is suitable also for **LEFT-HANDED** doctors. Adapting the unit to left-handed users is very quick and easy.

2 CONF IGUR ATION



CONF IGUR ATION

FIT YOUR SPACE

pick what you need

You can choose which kind of **INSTRUMENT TRAY** is more suitable for you, based on your habits. Also the choice of the **SPITTOON BOWL** is up to you: the dark glass spittoon bowl is removable to simplify the cleaning operation, The ceramic one is not removable.

HANGING TYPE



TOP MOUNTED



CERAMIC SPITTOON



GLASS SPITTOON

No matter what you choose,
there will be no additional costs.